

Figure 2: Melting-curve analysis for *ESR2*. Melting peaks of PCR products from leukocyte DNA and DNA treated with *Hpa*II, *Hha*I or *Sss*I methylases, resulting in methylation of 3, 7, 10, or 25 CpG sites.

methylase, or *Sss*I methylase (NEB) according to the manufacturer's recommendations.

For bisulfite conversion of DNA, the CpGenome DNA Modification Kit (Q-Biogene) was used, following the protocol provided by the manufacturer.

A 273-bp fragment of the estrogen receptor β gene (*ESR2*) promoter region corresponding to residues 1,912–2,185 (GenBank Accession No. AF191544) containing 25 CpGs was amplified in the LightCycler Instrument following bisulfite treatment (Figure 1). The PCR amplification was performed in a 10- μ l volume using the LightCycler FastStart DNA Master^{PLUS} SYBR Green I Kit, 20 ng bisulfite-treated DNA, 10 pmol primer Biseq-*ESR2*-fwd 5'-tt(ct)gtaggaggtagtgtgaag-3' and Biseq-*ESR2*-rev 5'-cacc(ag)accttaccactctaaaata-3'. The initial denaturing step at 95°C for 10 minutes was followed by 50 cycles of 95°C for 10 seconds, 58°C for 5 seconds and 72°C for 15 seconds. The melting-curve analysis was performed at 68°C according to the protocol supplied with the kit.

Results

We used a sequence from the estrogen receptor β gene (*ESR2*) CpG island to study the sensitivity of melting-curve analysis with the LightCycler Instrument. This *ESR2* sequence is hypermethylated in some prostate and breast cancers [5, 6]. As shown in Figure 1, the 273-base pair fragment amplified contains 25 CpG sites, three of which are located in *Hpa*II sites (CCGG) and seven are located in *Hha*I sites (GCGC). The sequence is completely unmethylated in leukocyte DNA from young healthy individuals.

DNA from healthy individuals was methylated *in vitro* using *Hpa*II methylase, *Hha*I methylase, a combination of these methylases, or *Sss*I methylase. These treatments introduce methylation at 3, 7, 10, or 25 sites. The melting-curve analysis of the PCR products from these DNAs is shown in Figure 2. Complete methylation shifted the melting point by 5.1°C from 83.7°C to 88.8°C. Partial methylation yielded less pronounced differences, but methylation of only three sites resulting in a shift of 0.7°C was well distinguished. Even partial methylation was sensitively detected.





Figure 3: Melting-curve analysis for *ESR2*. Melting peaks of PCR products from leukocyte DNA with different admixtures of DNA treated with *Hpa*II and *Hha*I methylases.



To ensure that melting-curve analysis with the LightCycler Instrument detected a fraction of methylated DNA against a background of unmethylated DNA, leukocyte DNA was mixed with various amounts of DNA methylated with *HpaII* and *HhaI* methylases (Figure 3). Admixture of 5–10% of methylated DNA resulted in detectable peak “shoulders” in the melting-curve analysis.

In contrast to DNA methylated by bacterial methylases *in vitro*, DNA from tissue or cell lines represents a mixture of DNAs with varying methylation patterns – if hypermethylation takes place at all. To assess the suitability of melting-curve analysis for such “real-world” samples, DNA from a range of cell lines was employed. These samples often yielded broader curves, as expected for DNA with varying degrees of hypermethylation (Figure 4). Some cell lines and tissues presented with a peak that clearly corresponds to unmethylated DNA (*e.g.*, the prostate cancer cell line PC3). Other peaks, however, extended towards the methylated control. For instance, the prostate cancer cell line LNCaP, known to contain a hypermethylated *ESR2* CpG island from bisulfite sequencing analysis, peaked at >86°C. Interestingly, and somewhat unexpectedly, several bladder cancer cell lines were found to harbor substantial degrees of hypermethylation, as exemplified by HT1376 and T-24 in Figure 4. According to this analysis, VmCub1 appears to contain unmethylated as well as methylated alleles.

Conclusion

The results indicate that melting-point analysis following bisulfite treatment and PCR using the LightCycler is well-suited to screening for DNA hypermethylation of CpG islands in cancer and other diseases. As few as three out of 25 sites and 5–10% methylated DNA were sufficient to detect hypermethylation in the exemplary sequence used. It may be possible to further increase the level of sensitivity by using shorter amplification products with higher

CpG content, different salt conditions, or other DNA stains. However, the present level of sensitivity should be quite sufficient to detect physiologically relevant hypermethylation in tumor tissues.

The high sensitivity of the technique is somewhat unexpected in view of theoretical considerations [4]. This is probably due to the fact that hypermethylation of DNA results in an increased GC content after bisulfite treatment. Therefore, PCR products from methylated DNA present as novel peaks at increased temperature compared with peaks obtained under normal circumstances. Any contaminant products, such as primer multimers, however, melt at lower temperatures. Therefore, the background level at higher temperatures is low, and even small shifts and minor shoulders can be relatively easily discerned.

In the present investigation, a simple screen revealed the presence of hypermethylated *ESR2* sequences in bladder cancer cell lines. Estrogen receptors are well-known to play a role in bladder function, but have hardly been studied in the context of bladder cancer. Thus, this screen raises the interesting question of whether such changes also occur in the respective tissues and may be relevant for tumorigenesis. ■

References

1. Jones PA, Baylin SB (2002) *Nat Rev Cancer* 3: 415–428
2. Clark SJ, Warnecke PM (2002) *Methods* 27: 99–100
3. Betz B *et al.* (2004) *Hum Mutat* 23: 612–620
4. Guldberg P, Worm J, Gronbaek K (2002) *Methods* 27: 121–127
5. Nojima D *et al.* (2001) *Cancer* 92: 2076–2083
6. Zhao C *et al.* (2003) *Oncogene* 22: 7600–7606

Product	Pack Size	Cat. No.
LightCycler 2.0 Instrument	1 instrument	03 531 414 201
FastStart DNA Master^{PLUS} SYBR Green I Kit	1 kit (96 reactions)	03 515 869 001

